



## PRODUCT BROCHURE



Affordable Lab Technology

# Contents

---

	<b>PAGE NO.</b>
Introduction	4
<b>T6U</b> UV-Vis Spectrophotometer	6
<b>T7</b> UV-Vis Spectrophotometer	10
<b>T7D</b> UV-Vis Spectrophotometer	14
<b>T8DCS</b> UV-Vis Spectrophotometer	18
<b>T9/T10</b> UV-Vis Spectrophotomete	22
<b>UVWIN</b> Software	24
<b>T3/T3M</b> Portable Spectrometer	28
Specifications	32
Accessories	35
Consumables	39
<b>A3</b> Atomic Absorption Spectrometer	40
<b>PF7</b> Atomic Fluorescence Spectrometer	44

Affordable Lab Technology



ACCURATE · RELIABLE · FLEXIBLE ·

**Since its development in the 1950's the UV-Visible Spectrophotometer has evolved into an accurate and reliable analytical tool and it has become one of the most utilised instruments in today's scientific laboratory.**

UV-Visible spectroscopy is a mature and well established analytical technique used extensively in many industry sectors including Environmental Analysis, Pharmaceutical Testing, Food and Beverage Production to name but a few. Persee Analytics manufacture an extensive range of UV and Visible Spectroscopy instrumentation guaranteed to meet the needs of your application. Further information on the UV-Vis product line along with a brief introduction to UV-Vis Spectroscopy can be found in this brochure.

#### **UV-Vis Spectroscopy**

UV-Vis Spectroscopy is an analytical method used to measure the absorbance of ultra-violet or visible radiation through an analyte. The molecular

absorption of the analyte corresponds to both excitation of valence electrons and excitation of electrons in different atomic orbitals.

UV-Vis Spectroscopy is an effective technique for both qualitative and quantitative analysis of organic and inorganic compounds.

UV-Vis Spectroscopy is based on the Lambert-Beer principle which states that the Absorbance of a solution (A) is directly proportional to its pathlength (l) and its concentration (c) when the wavelength of the incidence light remains fixed.

This is summarized in the following equation, where  $\epsilon$  is the molar absorptivity

$$A = \epsilon l c$$



## ACCESSIBLE · ANALYSIS

### UV-Vis Spectrophotometer

The UV-Vis Spectrophotometer is the analytical instrument used for the UV-Vis spectroscopic analysis. Spectrophotometers are available in different configurations however most can be categorized into either single beam, split beam or double beam types depending on the design of their optical system. Such types of instrument comprise the following components in their constructions:

- Light Source
- Monochromator
- Cell Compartment
- Detector
- Signal Processing System

### Split Beam Spectroscopy

The Split beam approach to UV-Vis Spectroscopy uses a single beam of light separated into Sample and

Reference beams by means of a beam splitter using a separate detector for each. Split beam instruments have a reference detector housed inside the instrument optics offering the advantage of optical stability as in double beam spectroscopy whilst using the single beam measurement technique.

### Double Beam Spectroscopy

The double beam approach to UV-Vis spectroscopy requires two beams of light, both having the same intensity to measure the Absorbance through sample and reference positions simultaneously. The Sample position is used for measurement of the analyte, whereas the reference position is used for the correction against a blank solution or sample matrix.

A clear advantage of the double beam optical system is the improvement in measurement stability and drift precision as a result of having a real-time feedback of both the reference and sample signals.

# T6

UV-VIS SPECTROPHOTOMETER



**The T6 is a high performance compact split beam spectrophotometer with a fixed 2nm spectral bandwidth.**

**The T6 range consists of two models:**

T6U (UV-Visible) operating within a wavelength range of 190 -1100nm. T6V (Visible) operating within a wavelength range of 325 -1100nm.

The instrument has a switched mode power supply accepting voltages in the range of 95 - 240V AC and supplied with either universal pathlength 5 cell changer or fixed path length 8 cell changer as standard.

The T6 delivers the functionality and accuracy of an advanced instrument at an affordable price.

**FEATURES & FUNCTIONS**

- High performance fixed 2nm spectral bandwidth.
- Low stray light 0.05%T (T60U).
- Wavelength accuracy +/- 1nm (T60U).
- Holographic blazed grating 1200lines/mm.
- Local control software for photometric fixed wavelength measurement.
- Easily upgraded to include quantitative analysis, multi wavelength spectrum & kinetics.
- Built in cell holder storage.
- Robust modular design with a small footprint.
- Can be used with UV-WIN software (optional).



Fixed pathlength 8-cell changer

# T6<sub>continued</sub>

## OPTICAL SYSTEM & COMPONENTS

High quality optical components ensure reliable analytical data with low stray light achieved using very low noise electronic circuits.

Deuterium and tungsten light sources deliver superior stability across the full wavelength range. Both types of lamp have self timers and are inexpensive and easy to replace when required.

The monochromator is completely sealed and the optical surfaces can be easily cleaned to maintain optimum reflectivity over the lifetime of the instrument.

A maintenance free high resolution direct stepper drive positions the grating precisely, which ensures reproducible wavelength scanning at different scan speeds, thus negating any wavelength peak shift.

The spectrophotometer shell is made from an environmentally friendly non corrosive material and a







# Specifications

T6U (UV-Visible)		T6V (Visible)
Optical system	Split beam ratio	Split beam ratio
Scan speed	Selectable	Selectable
Wavelength range	190 -1100nm	325 -1100nm
Wavelength accuracy	±1nm	± 2nm
Wavelength reproducibility	≤ 0.2nm	≤ 0.4nm
Spectral bandwidth	2nm	2nm
Photometric mode	Transmittance, Absorbance, Energy	Transmittance, Absorbance, Concentration
Photometric range	-0.3-3.0Abs	-0.3-3.0Abs
Photometric accuracy	0.002A (0 -0.5A), 0.004A (0.5 -1A) 0.3%T (0 -100%T)	0.002A (0 - 0.5A), 0.004A (0.5 - 1A) 0.3%T (0 -100%T)
Photometric reproducibility	0.001A (0 - 0.5A), 0.002A (0.5 - 1A) 0.15%T (0 -100%T)	0.001A (0 - 0.5A), 0.002A (0.5 - 1A) 0.15%T (0 -100%T)
Photometric noise	0.001A (500nm) 30min warm-up	0.001A (500nm) 30min warm-up
Baseline flatness	0.002A (200 -1000nm)	0.002A (325 -1000nm)
Baseline stability	0.001A/h (500nm,0Abs), 2hr warm-up	0.002A/h (500nm,0Abs), 2hr warm-up
Stray light	≤ 0.05%T (220nm NaI, 340nm NaNO <sub>2</sub> )	≤ 1.0%T (340nm NaNO <sub>2</sub> )
Standard Functionality	Photometric Measurement (Quantitative, Multi-wavelength, Spectrum and Kinetic measurements with program cards)	Photometric & Quantitative Measurement
Detector	Silicon photo diode	Silicon photo diode
Light source	Tungsten Halogen and Deuterium arc lamps	Tungsten Halogen
Display	Digital LCD display	Digital LCD Display
Printer	Mini Printer	Mini Printer
PC Interface	RS232	RS232
Software support	Local and UV-Win	Local and UV-Win
Power supply	Switch mode 95 - 250VAC 50 - 60Hz	Switch mode 95 - 250VAC 50 - 60Hz
Weight	11kg	11kg
Dimensions (Width, Depth, Height)	476(mm), 362(mm), 225(mm)	476(mm), 362(mm), 225(mm)

**Each Unit is supplied with the following as standard:**

- |   |                        |
|---|------------------------|
| 1 x Certificate of conformity                   | 1 x Power cord         |
| 1 Pair of quartz cells (T6U), glass cells (T6V) | 1 x Instruction manual |
| 1 x Quantitative program card (T6U Only)        | 1 x Dust cover         |
| 1 x Black block for dark current correction     | 1 x Packing list       |
| 1 x Fuse (2A)                                   |                        |

# T7

## UV-VIS SPECTROPHOTOMETER



660.0 nm 100.0 T

No.	T%	K * T%

Press START to measure

21:4



**The T7 is a high performance split beam spectrophotometer available with a fixed (2nm) or variable (0.5, 1, 2, 5nm) spectral bandwidth, which is innovative in terms of instrument application, mechanical and optical design, electronic control and software whilst retaining features that are well established and accepted through the industry.**

The T7 series of UV-Visible Spectrophotometers are able to carry out photometric measurement, spectrum scans, quantitative determination and DNA/Protein analysis. When interfaced to a PC using the UV-Win software, many more features are available including three dimensional spectrum, kinetic measurements, method and data storage, exportation of data in multiple formats and GLP administration features. Both instruments have a spectral range of 190-1100nm.

**The T7 range consists of two models:**

T7 UV-Vis instrument offering a fixed bandwidth of 2nm.

T74 UV-Vis instrument offering a variable bandwidth of 0.5, 1, 2 or 5nm.

**FEATURES & FUNCTIONS**

- High performance fixed (2nm) or variable (0.5, 1, 2, 5nm) spectral bandwidth.
- Wavelength accuracy +/- 0.3nm.
- Supplied with a motorised 8 cell changer and pre-aligned Tungsten and Deuterium lamps.
- Holographic blazed grating 1200 lines /mm.
- High degree of automation requiring minimal key depressions to start analysis.
- A number of optional accessories available which increase the flexibility of the analysis.
- Analysis for photometric measurement, spectrum scans, quantitative determination and DNA/Protein analysis.
- UV-Win software gives additional functionality including 3D spectrum analysis and compliance with GLP protocol.
- Simple mechanical structure and modular electronics make routine maintenance very easy.

# T7 continued




## OPTICAL SYSTEM & COMPONENTS

High quality optical components ensure reliable analytical data with low stray light achieved using very low noise electronic circuits.

The split beam ratio optics ensures good optical stability.

Pre-aligned Deuterium and tungsten light sources deliver superior stability across the full wavelength range. Both types of lamps are inexpensive and easy to replace when required.

The modular design allows easy access to all optical surfaces which can be easily cleaned to maintain optimum reflectivity over the lifetime of the instrument.



8 cell changer

# Specifications

Instrument Type	T7	T7S
Optical System	Split beam ratio	Split beam ratio
Scan Speed	Selectable	Selectable
Wavelength Range	190 -1100nm	190 -1100nm
Wavelength Accuracy	± 0.3nm	± 0.3nm
Wavelength Reproducibility	≤ 0.2nm	≤ 0.2nm
Spectral Bandwidth	2 nm	0.5, 1.0, 2.0, 5.0nm
Photometric Mode	Transmittance, Absorbance, Energy Concentration	Transmittance, Absorbance, Energy Concentration
Photometric Range	-0.3-3.0Abs	-0.3-3.0Abs
Photometric Accuracy	0.002A (0 - 0.5A), 0.004A (0.5 -1.0A), 0.3%T (0 - 100%T)	0.002A (0 - 0.5A), 0.004A (0.5 - 1.0A), 0.3%T (0 -100%T)
Photometric Reproducibility	0.001A (0 - 0.5A), 0.002A (0.5 - 1.0A), 0.15%T (0 -100%T)	0.001A (0 - 0.5A), 0.002A (0.5 - 1.0A), 0.15%T (0 -100%T)
Photometric Noise	0.001A (500nm) 30min warm-up	0.001A (500nm) 30min warm-up
Baseline Flatness	0.002A (200 - 1000nm)	0.002A (200 - 1000nm)
Baseline Stability	0.001A/h (500nm, 0Abs), 2hr warm-up	0.001A/h (500nm, 0Abs), 2hr warm-up
Stray light	≤ 0.12%T (220nm NaI, 340nm NaNO <sub>2</sub> )	≤ 0.12%T (220nm NaI, 340nm NaNO <sub>2</sub> )
Standard Functionality	Photometric, Quantitative, Spectrum and DNA measurements.	Photometric, Quantitative, Spectrum and DNA measurements.
Cell Holder	Automatic 8 Cell changer	Automatic 8 cell changer
Detector	Silicon photo diode	Silicon photo diode
Light Source	Tungsten Halogen and Deuterium arc lamps	Tungsten Halogen and Deuterium arc lamps
Display	Digital LCD display	Digital LCD display
Printer	Not available	Not available
PC Interface	RS232/USB	RS232/USB
Software Support	Local and UV Win	Local and UV Win
Power Supply	Switchable 120 - 230VAC 50 - 60Hz	Switchable 120 - 230VAC 50 - 60Hz
Weight	25 Kg	25Kg
Dimensions (Width, Depth, Height)	520mm, 420mm, 230mm	520mm, 420mm, 230mm

**Each Unit is supplied with the following as standard:**

- |   |                        |
|---|------------------------|
| 1 x Certificate of conformity                         | 1 x Power cord         |
| 1 x 8 position 10mm path length motorised cell holder | 1 x Instruction manual |
| 1 Pair 10mm Quartz cells                              | 1 x Dust cover         |
| 1 x Black block for dark current correction           | 1 x Packing list       |
| 1 x Fuse (2A)   |                        |

# T7D

UV-VIS SPECTROPHOTOMETER





**The T7D is a high performance double beam spectrophotometer available with a fixed (2nm) or variable (0.5, 1, 2, 5nm) spectral bandwidth, which is innovative in terms of instrument application, mechanical and optical design, electronic control and software whilst retaining features that are well established and accepted through the industry.**

The T7D series of UV-Visible Spectrophotometers are able to carry out photometric measurement, spectrum scans, quantitative determination and DNA/Protein analysis. When interfaced to a PC using the UV-Win software, many more features are available including three dimensional spectrum, kinetic measurements, method and data storage, exportation of data in multiple formats and GLP administration features. Both instruments have a spectral range of 190-1100nm.

#### **The T7D range consists of two models:**


T7D UV-Vis instrument offering a fixed bandwidth of 2nm.

T7DS UV-Vis instrument offering a variable bandwidth of 0.5, 1, 2 or 5nm.

#### **FEATURES & FUNCTIONS**

- High performance fixed (2nm) or variable (0.5, 1, 2, 5nm) spectral bandwidth.
- Wavelength accuracy +/- 0.3nm.
- Supplied with a motorised 8 cell changer and pre-aligned Tungsten and Deuterium lamps.
- Holographic blazed grating 1200 lines /mm.
- High degree of automation requiring minimal key depressions to start analysis.
- A number of optional accessories available which increase the flexibility of the instrument.
- Analysis for photometric measurement, spectrum scans, quantitative determination and DNA/Protein analysis.
- UV-Win software gives additional functionality including 3D spectrum analysis and compliance with GLP protocol.
- Simple mechanical structure and modular electronics make routine maintenance very easy.

# T7D continued



## OPTICAL SYSTEM & COMPONENTS

High quality optical components ensure reliable analytical data with low stray light achieved using very low noise electronic circuits.

The double beam optics ensure good optical stability. Pre-aligned Deuterium and Tungsten light sources deliver superior stability across the full wavelength range. Both types of lamps are inexpensive and easy to replace when required.

The modular design allows easy access to all optical surfaces which can be easily cleaned to maintain optimum reflectivity over the lifetime of the instrument.



# Specifications


Instrument Type	T7D	T7DS
Optical System	Double beam	Double beam
Scan Speed	Selectable	Selectable
Wavelength Range	190 -1100nm	190 -1100nm
Wavelength Accuracy	± 0.3nm	± 0.3nm
Wavelength Reproducibility	δ 0.2nm	δ 0.2nm
Spectral Bandwidth	2 nm	0.5, 1.0, 2.0, 5.0nm
Photometric Mode	Transmittance, Absorbance, Energy Concentration	Transmittance, Absorbance, Energy Concentration
Photometric Range	-0.3 - 3.0Abs	-0.3 - 3.0Abs
Photometric Accuracy	0.002A (0 - 0.5A), 0.004A (0.5 - 1.0A), 0.3%T (0 -100%T)	0.002A (0 - 0.5A), 0.004A (0.5 - 1.0A), 0.3%T (0 -100%T)
Photometric Reproducibility	0.001A (0 - 0.5A), 0.002A (0.5 - 1.0A), 0.15%T (0 -100%T)	0.001A (0 - 0.5A), 0.002A (0.5 -1.0A), 0.15%T (0 - 100%T)
Photometric Noise	0.001A (500nm) 30min warm-up	0.001A (500nm) 30min warm-up
Baseline Flatness	0.0015A (200 -1000nm)	0.0015A (200 -1000nm)
Baseline Stability	0.0008A/h (500nm, 0Abs), 2hr warm-up	0.0008A/h (500nm, 0Abs), 2hr warm-up
Stray light	≤ 0.12%T (220nm NaI, 340nm NaNO <sub>2</sub> )	≤ 0.12%T (220nm NaI, 340nm NaNO <sub>2</sub> )
Standard Functionality	Photometric, Quantitative, Spectrum and DNA measurements.	Photometric, Quantitative, Spectrum and DNA measurements.
Cell Holder	Automatic 8 Cell changer	Automatic 8 cell changer
Detector	Silicon photo diode	Silicon photo diode
Light Source	Tungsten Halogen and Deuterium arc lamps	Tungsten Halogen and Deuterium arc lamps
Display	Digital LCD display	Digital LCD display
Printer	Not available	Not available
PC Interface	RS232/USB	RS232/USB
Software Support	Local and UV Win	Local and UV Win
Power Supply	Switchable 120 - 230VAC 50 - 60Hz	Switchable 120 - 230VAC 50 - 60Hz
Weight	25 Kg	25 Kg
Dimensions (Width, Depth, Height)	520mm, 420mm, 230mm	520mm, 420mm, 230mm

**Each Unit is supplied with the following as standard:**

- |   |                        |
|---|------------------------|
| 1 x Certificate of conformity                         | 1 x Fuse (2A)          |
| 1 x 8 position 10mm path length motorised cell holder | 1 x Power cord         |
| 1 Pair 10mm Quartz cells                              | 1 x Instruction manual |
| 1 x Black block for dark current correction           | 1 x Packing list       |

# T8DCS

UV-VIS SPECTROPHOTOMETER



**The T8DCS is a high performance double beam spectrophotometer with a variable spectral bandwidth from 0.1- 5nm, selected by a continuous variable slit.**

The Czerny-Turner monochromator with a holographic grating keeps stray light to a minimum and offers excellent optical resolution. The use of a photomultiplier tube as a detector offers exceptional sensitivity.

The T8DCS' true double beam optical system coupled with an efficient and well proven electronic control system ensures high stability and low background noise.

#### **FEATURES & FUNCTIONS**

- Photomultiplier tube detection provides exceptional sensitivity.
- Wavelength accuracy  $\pm 0.3\text{nm}$  (Automatic Wavelength Correction).
- User selectable spectral bandwidth between 0.1-5nm.
- User friendly design allows easy light source replacement and routine maintenance.
- Sample compartment design enables use of a wide range of optional accessories.
- UV-WIN software offers many operational and data processing capabilities and is supplied as standard with the T8DCS.

# T8DCS<sub>continued</sub>



## OPTICAL SYSTEM & COMPONENTS

The T8DCS features an advanced continuous variable bandwidth feature making it the instrument of choice for applications with a demand for precise and accurate control of wavelength resolution. This feature allows the user to specify exactly what bandpass is required in the range of 0.1-5nm.

The double beam optical design combined with a high specification holographic grating gives excellent wavelength separation allowing the user to measure close adjacent wavelengths with excellent sensitivity.

The modular design of the sample compartment allows for ease of use of a wide range of optional accessories ensuring accurate analysis of various sample types including liquids, thin films and powders.

The user friendly design of the lamp compartment allows easy replacement and simplified routine maintenance of the Deuterium and Tungsten lamps.


# Specifications

Specifications	T8DCS
Optical System	Double beam
Scan Speed	Selectable
Wavelength Range	190 - 900nm
Wavelength Accuracy	± 0.3nm
Wavelength Reproducibility	≤ 0.1nm
Spectral Bandwidth	Continuous slit 0.1 - 5.0nm with 0.1nm interval
Photometric Mode	Transmittance, Absorbance, Energy Concentration, All Using UVWin Software
Photometric Range	-4.0 - 4.0Abs
Photometric Accuracy	0.002A (0 - 0.5A), 0.004A (0.5 - 1.0A), 0.3%T (0 -100%T)
Photometric Reproducibility	0.001A (0 - 0.5A)
Photometric Noise	0.0004A (500nm) 30min warm-up
Baseline Flatness	0.001A (200 - 850nm)
Baseline Stability	0.0008A/h (500nm, 0Abs), 2hr warm-up
Stray light	≤ 0.01 %T (220nm NaI, 340nm NaNO <sub>2</sub> )
Standard Functionality	No stand alone function
Cell Holder	Fixed position sample and reference
Detector	Photo multiplier tube
Light Source	Tungsten Halogen and Deuterium arc lamps
Display	No display
Printer	Not available
PC Interface	RS232/USB
Software Support	UV Win
Power Supply	Switchable 120 - 230VAC 50 - 60Hz
Weight	43Kg
Dimensions (Width, Depth, Height)	545mm, 580mm, 270mm

**Each Unit is supplied with the following as standard:**

- |   |                        |
|---|------------------------|
| 1 x Certificate of conformity                                     | 1 x Fuse (2A)          |
| 1 x Standard fixed position cell holder<br>(sample and reference) | 1 x Power cord         |
| 1 x Pair Quartz cells   | 1 x Instruction manual |
| 1 x Black block for dark current correction                       | 1 x Dust cover         |
| 1 x UV Win Software disk  | 1 x Packing list       |

# T9DCS/T10DCS



**The T9DCS/T10DCS Series Spectrophotometer incorporates dual monochromator technology making it well suited to even the most demanding of applications in all areas of UV-Visible Spectroscopy including:**

- Pharmaceutical
- Metrological Verification
- Food Safety
- Material Science
- Biotechnology
- Research



Wireless instrument control and data acquisition can be achieved from a Tablet PC using Wi-Fi technology allowing the analyst to move freely around the laboratory whilst also keeping track of sample measurements.

The optical design of both the T9DCS and T10DCS offer extremely low stray light characteristics ( $\leq 0.00004\%T$  NaI, 220 nm) which allows for an extensive photometric range (-8.0 – 8.0Abs). Measurements at deep ultra-violet wavelengths can also be achieved with use of Nitrogen purged optics.

The instrument can be optically configured to suit the needs of the sample by means of a continually adjustable slit for precise control of spectral resolution and beam size adjustment by means of an attenuating wheel.

Precise wavelength accuracy is ensured by the integrated Mercury Emissions Lamp used for automatic correction of spectral deviation.

A whole host of specialised accessories are available to suit the specific requirement of the sample, these include:

- Both 60mm and 150mm Integrating Sphere for Diffuse reflectance measurements.
- Absolute, and Specular reflectance measurements accessories.
- Polarizing Optics.
- Thermostatic Cell Holders for temperature control.
- Various long and short pathlength cell holders.
- Automated cell changers for both sample and reference beams.
- Tablet dissolution accessory for pharmaceutical quality control.

## Specifications Further product information available soon.

Specifications	T9DCS	T10DCS
Optical System	Dual Monochromator Double Beam	Dual Monochromator Double Beam
Light Source	D2 Lamp – UV Region W Lamp – Visible Region Hg Lamp – Wavelength Correction	D2 Lamp – UV Region W Lamp – Visible Region Hg Lamp – Wavelength Correction
Wavelength Range	185–900nm	185–900nm
Wavelength Accuracy	$\pm 0.2\text{nm}$	$\pm 0.2\text{nm}$
Wavelength Reproducibility	$\leq 0.1\text{nm}$ (D2 lamp)	$\leq 0.1\text{nm}$ (D2 lamp)
Spectral Bandwidth	0.1 – 5nm Continually Adjustable	0.1 – 5nm Continually Adjustable
Stray Light	$\leq 0.0001\%T$ (NaI, 220 nm) $\leq 0.0001\%T$ (NaNO <sub>2</sub> , 360 nm)	$\leq 0.00004\%T$ (NaI, 220 nm) $\leq 0.00002\%T$ (NaNO <sub>2</sub> , 360 nm)
Photometric Range	-6.0Abs–6.0Abs	-8.0Abs–8.0Abs
Photometric Accuracy	$\pm 0.004A$ @2.0A $\pm 0.003A$ @1.0A $\pm 0.002A$ @0.5A $\pm 0.3\%$	$\pm 0.004A$ @2.0A $\pm 0.003A$ @1.0A $\pm 0.002A$ @0.5A $\pm 0.3\%$
Photometric Reproducibility	$\leq 0.002A$ @2.0A $\leq 0.0008A$ @1.0A $\leq 0.0004A$ @0.5A $\leq 0.1\%$	$\leq 0.002A$ @2.0A $\leq 0.0008A$ @1.0A $\leq 0.0004A$ @0.5A $\leq 0.1\%$
Baseline Flatness	$\pm 0.0008\text{Abs}$	$\pm 0.0005\text{Abs}$
Noise	0% Noise: $\leq 0.01\%$ ; 100% T Noise: $\leq 0.1\%$ ;	0% Noise: $\leq 0.01\%$ ; 100% T Noise: $\leq 0.1\%$ ;
Communication port	RS232C, USB, Wifi	RS232C, USB, Wifi

We reserve the right to modify, revise/upgrade, suspend or discontinue any Product in whole or in part, either temporarily or permanently, with or without notice.

# UVWin



**UVWin is a powerful, intuitive Software product used for connectivity to the range of bench top UV-Vis Spectrophotometers.**

The UVWin software offers complete instrument control along with data acquisition and a whole host of mathematical tools for interpretation of measurement results. The UVWin software is separated into four key workspaces:

- Spectral Analysis
- Quantitative Analysis
- Kinetic Analysis
- Photometric Analysis

## SPECTRUM WORKSPACE


- Use the spectrum workspace to scan across a user-defined spectral range measuring in either absorbance or transmission.
- Use the "Peak Pick" tool to determine the wave-length at which peaks and valleys have occurred whilst also being able to determine their amplitude.
- View spectral overlay in the 3D display mode.
- Perform 1st, 2nd, 3rd and 4th order differentiation on sample scans for Derivative Spectroscopy.
- Export measurement data into Word, Excel, CSV and ASCII formats.
- Create method files for routine analysis whilst also being able to save measurement data.






## QUANTITATIVE WORKSPACE

- Use the Quantitative workspace to determine the concentration of unknown samples.
- Create a calibration curve using a series of standard solution or by directly entering the coefficients for the calibration curve equation.
- Specify 1st, 2nd, 3rd and 4th order correlation for the best calibration curve fit.
- Set Quality Control monitors to take user specified action in the event of measurement results falling outside user defined limits.
- Export measurement data into Word, Excel, CSV and ASCII formats.
- Create method files for routine analysis whilst also being able to save measurement data.



## KINETIC WORKSPACE

- Monitor the change of Absorbance or Transmission as a function of time for Enzyme type reactions.
- Use in conjunction with a Flowcell for sample introduction or Peltier water circulator for temperature control.
- Specify data intervals and acquisition time for up to 24 hour reactions.
- Export measurement data into Word, Excel, CSV and ASCII formats.
- Create method files for routine analysis whilst also being able to save measurement data.










# UVWinGLP


UVWin GLP offers all of the features and functionality of UVWin whilst also offering extensive Administrative capabilities along with a detailed audit trail.

## ADMINISTRATION

- Administrative settings can be made where Analysts may require conformity to GLP/GMP/GRP.
- Create User groups specifying exactly what actions they are able to perform.
- Add New Users to custom User Groups to determine their privilege settings.
- Automatically log software activity in an Audit Trail.
- Use Password control to ensure Users are logged in for instrument usage.



Date	Time	User	Action
2011-01-25	09:47:15	Administrator	Login
2011-01-25	09:48:12	Administrator	Open Data File to C:\Program Files\U...
2011-01-25	09:48:4	Administrator	Open Data File to C:\Program Files\U...
2011-01-25	09:49:05	Administrator	Open Data File to C:\Program Files\U...
2011-01-25	09:50:3	Administrator	Login
2011-01-25	09:50:55	Administrator	Open Data File to C:\Program Files\U...
2011-01-25	10:01:2	Administrator	Open Data File to C:\Program Files\U...
2011-01-25	10:02:0	Administrator	Open Data File to C:\Program Files\U...
2011-01-25	10:04:2	Administrator	Open Data File to C:\Program Files\U...
2011-01-25	10:05:4	Administrator	View Log



## CERTIFICATION

UVWin GLP has been evaluated and tested by a third party software validation specialist. As a result it was found that UVWin GLP offers all of the features and functions required for use in compliance with the guidance specified in:

- 21CFR Part 11- Electronic Records; Electronic Signatures
- Guidance for Industry Part 11, Electronic Records; Electronic Signatures — Scope and Application, August 2003

# T3/T3M continued

## PORTABLE SPECTROMETER



The T3 Portable Spectrometer is a compact portable instrument based on advanced CCD detection technology. Whilst being highly compact the instrument boasts all of the features of a conventional bench-top Spectrophotometer including Spectrum Scanning, Photometric, Quantitative and Kinetic Methods.

A Windows CE embedded operating system and touch screen TFT interface allow for ease of use and extensive data storage. Field measurement data can be transferred from internal instrument storage to a PC via USB connection and the T3 Data Viewer Software which offers an extensive tool set for data interpretation and reporting. The T3 Spectrometer is supplied with a rugged carry case and a wide range of accessories suited to your sampling requirements.

The T3M differentiates itself from the standard T3 by means of pre-programmed test methods for the Spectroquant® series of reagent test kits from Merck Chemicals.

Persee Analytics have selected the Spectroquant® Test Kits from Merck Chemicals as a partner in order to offer an application specific instrument, targeted specifically for environmental measurements. The T3M Spectrometer is supplied pre-programmed with calibration data for each of the Spectroquant® test methods.




Spectroquant® test kits can offer an analytical solution for the following parameters by means of 130 different test kits:

- Drinking water
- Surface water
- Process water
- Municipal or industrial wastewater
- Beverages
- Disinfectant control


Note: Spectroquant is a registered trademark of Merck Chemicals, Germany.



Windows CE Based Operating System with 320 x 240 touch screen TFT Interface.



The Spectrum workspace allows high speed spectral scanning, with zoom and peak identification tools. Spectral Scans can be performed in the field, stored to instrument memory and later transferred to the




curves, and measure concentration of unknown samples. Curves can be constructed in 1st – 4th order, whilst both methods and measurement data can be saved to instrument memory.

Setting	Photometry (unsaved)				Brows
Sel.	ID	Mode	500.0nm	K <sup>+</sup>	
	1	Abs	22.000	44.000	
	2	Abs	51.000	102.000	
	3	Abs	36.000	72.000	
	4	Abs	10.000	20.000	
	5	Abs	72.000	144.000	

Use the Photometric workspace to quickly and easily

Sel	Name	Conc.	Abs	Measure Time	Sample Info
	sample 1	2.93mg/L	0.117	03/02/2012 10:11	;
	sample 2	2.31mg/L	0.409	10/02/2012 10:58	;
	sample 3	2.39mg/L	0.372	10/02/2012 10:58	;




Absorbance or Transmission. Set a K factor where multiplications are required to determine sample concentration. Once the measurement is complete store to instrument memory for future recall. The Kinetic workspace enables the measurement of Absorbance or Transmission as a function of time. Use the zoom and peak pick features to obtain a better view of the Kinetic curve. Measurement data can be saved and recalled at any time.





# T3/T3M



User and admin rights are easily controlled from the GLP feature in the settings menu. Create user groups and specify their privilege level then add new users to a specified group.




Use the universal cell holder to measure various pathlength rectangular cells and rounds test tubes accommodating all of your sampling requirements.



Control instrument configuration from the settings menu.



Use the fibre dip probe for in-situ sample measurements.



Once all of the required field analysis has been performed and measurement data stored to instrument memory the Spectrometer can be connected to the s3 Data Viewer Software via USB for transfer of analysis data from all of the instrument workspaces. Use the Data Viewer Software to further interpret analysis results, export data into a wide variety of formats and produce analysis reports for storage or printing.

# Specifications

<b>Light Source</b>	Convergent tungsten lamp with 7000hr lifespan.	
<b>Measurement Workspaces</b>	Spectrum Quantitative Kinetics Photometric	
<b>Optical System</b>	Polychromatic with concave holographic grating	
<b>Detector</b>	CCD Sony ILX511 2048 Pixels	
<b>Sampling Accessories</b>	Fibre Dip Probe with 10mm and 20mm pathlength tips (optical accessory) Rectangular Cuvette Holder Cylindrical Test Tube Holder	
<b>Power supply</b>	Built-in re-chargeable battery with 5 hr usage	
<b>Printer</b>	Micro printer (optional)	
<b>Operating system</b>	Windows Embedded CE 6.0 with 2GB Flash Memory	
<b>Input/display</b>	320 x 240 True Colour TFT Touch Screen	
<b>Specifications</b>	Wavelength range	380nm – 800nm
	Wavelength resolution	0.4nm
	Spectral bandwidth	4±0.8nm
	Wavelength accuracy	±1.0 nm
	Wavelength repeatability	~ 0.1nm
	Baseline flatness	±0.005Abs
	Noise	0.5%
	Drift	1.0%
	Stray light	0.5%
	Photometric accuracy	±1.0%
	Photometric repeatability	0.3%
<b>Dimension</b>	280 x 170 x 110mm	
<b>Environmental temperature</b>	Operating environment 5-30°C Storage environment -20-55°C	

# Specifications

	T6U (UV-Visible)	T6V (Visible)	T7
Optical System	Split beam ratio	Split beam ratio	Split beam ratio
Scan Speed	Selectable	Selectable	Selectable
Wavelength Range	190 -1100nm	325 -1100nm	190 -1100nm
Wavelength Accuracy	± 1nm	± 2nm	± 0.3nm
Wavelength Reproducibility	< 0.2nm	< 0.4nm	< 0.2nm
Spectral Bandwidth	2nm	2nm	2nm
Photometric Mode	Transmittance, Absorbance, Energy	Transmittance, Absorbance, Concentration	Transmittance, Absorbance, Energy and Concentration
Photometric Range	-0.3 - 3.0Abs	-0.3 - 3.0Abs	-0.3 - 3.0Abs
Photometric Accuracy	0.002A (0 - 0.5A), 0.004A (0.5 -1.0A), 0.3%T (0 -100%T)	0.002A (0 - 0.5A), 0.004A (0.5 -1.0A), 0.3%T (0 -100%T)	0.002A (0 - 0.5A), 0.004A (0.5 -1.0A), 0.3%T (0 -100%T)
Photometric Reproducibility	0.001A (0 - 0.5A), 0.002A (0.5 - 1.0A), 0.15%T (0 -100%T)	0.001A (0 - 0.5A), 0.002A (0.5 -1.0A), 0.15%T (0 -100%T)	0.001A (0 - 0.5A), 0.002A (0.5 - 1.0A), 0.15%T (0 -100%T)
Photometric Noise	0.001A (500nm) 30min warm-up	0.002A (500nm) 30min warm-up	0.001A (500nm) 30min warm-up
Baseline Flatness	0.002A (200 -1000nm)	0.002A (325 -1000nm)	0.002A (200 -1000nm)
Baseline Stability	0.001A/h (500nm, 0Abs), 2hr warm-up	0.002A (500nm, 0Abs), 2hr warm-up	0.001A/h (500nm, 0Abs), 2hr warm-up
Stray light	< 0.05%T (220nm NaI, 340nm NaNO <sub>2</sub> )	< 1.0%T (340nm NaNO <sub>2</sub> )	< 0.12%T (220nm NaI, 340nm NaNO <sub>2</sub> )
Integration Time	N/A	N/A	N/A
Standard Functionality	Photometric Measurement. (Quantitative, Multi-wavelength, Spectrum and Kinetic measurements with	Photometric and Quantitative measurements.	Photometric, Quantitative, Spectrum and DNA measurements.
Cell Holder	Automatic 8 Cell Changer	Automatic 8 Cell Changer	Automatic 8 Cell changer
Detector	Silicon photo diode	Silicon photo diode	Silicon photo diode
Light Source	Tungsten Halogen and Deuterium arc lamps	Tungsten Halogen lamps	Tungsten Halogen and Deuterium arc lamps
Display	Digital LCD display	Digital LCD display	Digital LCD display
Printer	Mini printer	Mini printer	N/A
PC Interface	RS232	RS232	RS232/USB
Software Support	Local and UV Win	Local and UV Win	Local and UV Win
Power Supply	Switch mode 95 - 250VAC 50 - 60Hz	Switch mode 95 - 250VAC 50 - 60Hz	Switchable 120 – 230VAC 5 60Hz
Weight	11Kg	11Kg	25Kg
Dimensions (Width, Depth, Height)	476mm, 362mm, 225mm	476mm, 362mm, 225mm	520mm, 420mm, 230mm

T7S	T7D	T7DS	T8DCS
Split beam ratio	Double beam	Double beam	Double beam
Selectable	Selectable	Selectable	Selectable
190 -1100nm	190 -1100nm	190 -1100nm	190 - 900nm
± 0.3nm	± 0.3nm	± 0.3nm	± 0.3nm
< 0.2nm	< 0.2nm	< 0.2nm	< 0.1nm
0.5, 1.0, 2.0, 5.0nm	2nm	0.5, 1.0, 2.0, 5.0nm	Continuous slit 0.1 - 5.0nm (0.1nm interval)
Transmittance, Absorbance, Energy and Concentration	Transmittance, Absorbance, Energy and Concentration	Transmittance, Absorbance, Energy and Concentration	Transmittance, Absorbance, Energy and Concentration
-0.3 - 3.0Abs	-0.3 - 3.0Abs	-0.3 - 3.0Abs	-4.0 - 4.0Abs
0.002A (0 - 0.5A), 0.004A (0.5 -1.0A), 0.3%T (0 -100%T)	0.002A (0 - 0.5A), 0.004A (0.5 - 1.0A), 0.3%T (0 -100%T)	0.002A (0 - 0.5A), 0.004A (0.5 -1.0A), 0.3%T (0 -100%T)	0.002A (0 - 0.5A), 0.004A (0.5 - 1.0A), 0.3%T (0 - 100%T)
0.001A (0 - 0.5A), 0.002A (0.5 - 1.0A), 0.15%T (0 -100%T)	0.001A (0 - 0.5A), 0.002A (0.5 -1.0A), 0.15%T (0 -100%T)	0.001A (0 - 0.5A), 0.002A (0.5 -1.0A), 0.15%T (0 -100%T)	0.001A (0 - 0.5A)
0.001A (500nm) 30min warm-up	0.001A (500nm) 30min warm-up	0.001A (500nm) 30min warm-up	0.0004A (500nm) 30min warm-up
0.002A (200 – 1000nm)	0.0015A (200 -1000nm)	0.0015A (200-1000nm)	0.001A (200 - 850nm)
0.001A/h (500nm, 0Abs), 2hr warm-up	0.0008A/h (500nm, 0Abs), 2hr warm-up	0.0008A/h (500nm, 0Abs), 2hr warm-up	0.0008A/h (500nm, 0Abs), 2hr warm-up
< 0.12%T (220nm NaI, 340nm NaNO <sub>2</sub> )	< 0.12%T (220nm NaI, 340nm NaNO <sub>2</sub> )	< 0.12%T (220nm NaI, 340nm NaNO <sub>2</sub> )	< 0.01%T (220nm NaI, 340nm NaNO <sub>2</sub> )
N/A	N/A	N/A	N/A
Photometric, Quantitative, Spectrum and DNA measurements.	Photometric, Quantitative, Spectrum and DNA measurements.	Photometric, Quantitative, Spectrum and DNA measurements.	No stand alone function
Automatic 8 cell changer	Automatic 8 Cell changer	Automatic 8 cell changer	Fixed position sample and reference
Silicon photo diode	Silicon photo diode	Silicon photo diode	Photo multiplier tube
Tungsten Halogen and Deuterium arc lamps	Tungsten Halogen and Deuterium arc lamps	Tungsten Halogen and Deuterium arc lamps	Tungsten Halogen and Deuterium arc lamps
Digital LCD display	Digital LCD display	Digital LCD display	No display
N/A	N/A	N/A	N/A
RS232/USB	RS232/USB	RS232/USB	RS232/USB
Local and UV Win	Local and UV Win	Local and UV Win	UV Win
Switchable 120 – 230VAC 50 – 60Hz	Switchable 120 - 230VAC 50 - 60Hz	Switchable 120 - 230VAC 50 - 60Hz	Switchable 120 - 230VAC 50 - 60Hz
25Kg	25Kg	25Kg	43Kg
520mm, 420mm, 230mm	520mm, 420mm, 230mm	520mm, 420mm, 230mm	545mm, 580mm, 270mm

# Specifications

Model	T9DCS	T10DCS	T3/T3M
Optical System	Double beam	Double beam	Polychromatic with concave holographic grating
Scan Speed	1000nm/min 1500nm/min location speed	1000nm/min 1500nm/min location	-
Wavelength Range	185-900nm	185-900nm	380-800nm
Wavelength Accuracy	±0.2nm	±0.2nm	± 1.0nm
Spectral Bandwidth	±20%(0.1nm-5nm continuously selectable)	±20%(0.1nm-5nm continuously selectable)	4nm ± 0.8nm
Mode	Transmittance, Absorbance, Energy, Concentration, Absorbance	Transmittance, Concentration	Transmittance, Energy and Concentration
Photometric Range	-6.0Abs-6.0Abs	-8.0Abs-8.0Abs	-
Photometric Accuracy	±0.004A @2.0 A ±0.003A @1.0 A ±0.002A @0.5 A ±0.3%	±0.004A @2.0 A ±0.003A @1.0 A ±0.002A @0.5 A ±0.3%	± 1.0%
Photometric Reproducibility	<0.002A @2.0 A <0.0008A / 1.0 A <0.0004A/ 0.5 A <0.1%	<0.002A @2.0 A <0.0008A / 1.0 A <0.0004A/ 0.5 A <0.1%	≤ 0.3%
Photometric Noise	0% Noise: <0.01% 100% Noise: <0.1%	0% Noise: <0.01% 100% Noise: <0.1%	-
Baseline Flatness	±0.0008Abs	±0.0005Abs	± 0.005Abs
Baseline Stability	0.1%/hr	0.1%/hr	-
Stray light	<0.0001%T (NaI, 220 nm) <0.0001%T (NaNO <sub>3</sub> , 360 nm)	<0.00004%T (NaI, 220 nm) <0.00002%T (NaNO <sub>3</sub> , 360 nm)	≤ 0.5%
Integration Time	N/A	N/A	-
Standard Functionality			Spectrum, Quantitative, Kinetics, Photometric
Cell Holder	Automatic 8 Cell Changer	Automatic 8 Cell Changer	-
Detector	Photo multiplier tube	Photo multiplier tube	CCD Sony ILX511 2048 Pixels
Light Source	D2 lamp, Tungsten lamp, Mercury lamp for WL correction.	D2 lamp and Tungsten lamp. Mercury lamp for WL correction.	Convergent tungsten lamp with 7000hr lifespan
Display	N/A	N/A	320 x 240 True Colour TFT Touch Screen
Printer	N/A	N/A	Micro printer (optional)
Pc Interface	RS232C,USB , Wifi	RS232C,USB ,Wifi	-
Software Support	UV Win	UV Win	-
Power Supply	Switchable 120-230VAC 50-60Hz	Switchable 120-230VAC 50-60Hz	Built-in re-chargeable battery with 5 hr usage
Weight	62Kg	62Kg	-
Dimensions (width, depth, height)	913mm, 645mm, 298mm	913mm, 645mm, 298mm	280mm, 170mm, 110mm

# Accessories

I 6



Available Program cards include:

- Quantitative card 21602-2801-00
- Spectrum scanning/Kinetics card 21604-2801-00
- Multi wavelength card 21605-2801-00
- Palm Oil Analysis Card (DOBI)



- CH16-1**  
**Constant temperature cell holder**
- Cell Pathlength: 10mm
  - Number of Cells: 5
  - Requires PTC-2 Peltier Water Circulator



- PS16-2**  
**Sipper Pump**
- Pump Speed: 0.1 - 250 RPM
  - Speed resolution: 0.1 RPM less than 30 RPM Speed and 1 RPM above 30 RPM Speed



- DS16-1**  
**Angle adjustable sample holder**
- Maximum angle: 45 Degrees
  - Minimum Sample Size: 4mm (Width)
  - Maximum Sample Size: 80 x 55 x 5mm



- PTC-2**  
**Peltier**
- Temperature range: 5 - 75°C
  - Use in conjunction with CH16-1



- TR16-1**  
**Test tube holder**
- Test Tube Diameter: 15 - 25mm
  - Test Tube Height: 90 - 120mm



- USB Printer driver P2U**
- Connect the T60 to specific USB printers

T 7



- LS181-1**  
**5 cell holder**
- Cell Pathlength: 5 - 50mm (adjustable)
  - Number of Cells: 5



- CH181-1**  
**Constant temperature sample holder**
- Cell Pathlength: 10mm
  - Number of Cells: 5
  - Requires PTC-2 Peltier Water Circulator



- DS181-1**  
**Angle adjustable holder**
- Maximum angle: 45 Degrees
  - Minimum Sample Size: 4mm (Width)
  - Maximum Sample Size: 80x55x5mm



- TR181-1**  
**Test tube holder**
- Test Tube Diameter: 15 - 25mm
  - Test Tube Height: 90 -120mm



- S181-1**  
**Solid Sample Holder**
- Maximum sample size: 80mm x 55mm x 5mm

**8F** offer a complete range of cuvettes

## T7 continued



- MR181-1**  
**Specular Reflection**
- Incidence angle: 5°
  - Size of Sample Area Measured: 11 × 9mm to 60 × 40mm
  - Spectral Range: 200 - 1100nm



- PS181-2**  
**Sipper Pump**
- Pump Speed: 0.1 - 250 RPM
  - Speed resolution: 0.1 RPM less than 30 RPM Speed and 1 RPM above 30 RPM Speed



- MH181-1**  
**Micro cell holder**
- Pathlength: 10mm
  - Minimum Cell Window Width: 2mm
  - Minimum Cell Window Height: 10mm



- PTC-2**  
**Peltier**
- Temperature range: 5 - 75°C
  - Use in conjunction with CH181-1

## T7D



- LS188-1**  
**5 cell holder**
- Cell Pathlength: 5-50mm (adjustable)
  - Number of Cells: 5



- CH188-1**  
**Constant temperature holder**
- Cell Pathlength: 10mm
  - Number of Cells: 2 Cells (one for Sample and one for Reference)
  - Requires PTC-2 Peltier Water Circulator



- PS181-2**  
**Sipper Pump**
- Pump Speed: 0.1 - 250 RPM
  - Speed resolution: 0.1 RPM less than 30 RPM Speed and 1 RPM above 30 RPM Speed



- PTC-2**  
**Peltier**
- Temperature range: 5 - 75°C
  - Use in conjunction with CH188-1

## T8DCS



- CH92**  
**Constant temperature holder**
- Cell Pathlength: 10mm
  - Number of Cells: 2 Cells (one for Sample and one for Reference)
  - Requires PTC-2 Peltier Water Circulator



- S19-1**  
**Solid Sample Holder**
- Maximum Sample Size: 80mm × 55mm × 5mm
  - Sample and reference beams



- IS92**  
**Integrating Sphere**
- Incidence angle: Sample 0° reference 8°
  - Minimum Sample Size for Diffuse Reflectance: 15mm × 25mm
  - Minimum Sample Size for Transmission: 20mm Diameter
  - Wavelength Range: 230-850nm with a 5nm Bandpass
  - Sphere Diameter: 58mm



- DS19-1**  
**Angle adjustable cell holder**
- Maximum angle: 45 Degrees
  - Minimum Sample Size: 4mm (Width)
  - Maximum Sample Size: 80 × 55 × 5mm

**8F offer a complete range of cuvettes**

## T8%\$4 continued




- TR19-1**  
**Test tube holder**
- Test Tube Diameter: 15 - 25mm
  - Test Tube Height: 90 - 120mm




- MR19-1**  
**Specular reflection**
- Incidence angle: 5°
  - Size of Sample Area Measured: 30 x 15mm
  - Spectral Range: 200 - 900nm



- LS19-1**  
**Long path-length holder**
- Pathlength: 5, 10, 20, 30, 40, 50, 100mm
  - Number of Cells: 2 Cells (one for Sample and one for Reference)



- MH19-1**  
**Micro cell holder**
- Pathlength: 10mm
  - Minimum Cell Window Width: 2mm
  - Minimum Cell Window Height: 10mm



- MH19-2**  
**Ultra micro cell holder**
- Pathlength: 10mm
  - Minimum Cell Window Width: 2mm
  - Minimum Cell Window Height: 5mm



- PS92**  
**Sipper Pump**
- Pump Speed: 0.1 - 250 RPM
  - Speed resolution: 0.1 RPM less than 30 RPM Speed and 1 RPM above 30 RPM Speed



- PTC-2**  
**Peltier**
- Temperature range: 5 - 75°C
  - Use in conjunction with CH19-1



- AS91**  
**8 Cell Changer**
- 8 position cell changer
  - Pathlength: 10mm


## T9DCS/T10DCS



Double beam 60mm integration sphere



Double beam 150mm integration sphere




Double beam relatively specular reflection accessories



Double beam thermostatic cell holder



Double beam test tube holder



Double beam micro cell holder



Adjustable solid sample holder




Double beam sipper pump accessories

**8F** offer a complete range of cuvettes



## T9DCS/T10DCS continued



Double beam 5-cell thermostatic holder



Double beam 50mm long path length 5-cell holder



Double beam 100mm long path length cell holder



Double beam dissolution accessories



Double beam 8-cell changer

## I3/I3M



**30-010  
Fibre Dip Probe**  
● 10mm and 20mm pathlength tips



Crushing module



Measurement module



Ultrasonic module



Temperature control module



Distillation module



Vortex module



Food Sample Extraction module



Centrifuge module



Arsenic hydride measurement module



Hydrogen sulfide measurement module



Titration module

**We offer a complete range of cuvettes**

## T3113M continued



Honey moisture analysis module



Module Package



Conductivity measurement module



Water quality analysis module



Soil sampling module



Universal module for soil and fertilizer analysis



Available boron analysis module



Basic module



Stirring module



Filtration module



Soil Sample Extraction module



Heating extraction module



-16 test tube module



-25 test tube module



Constant volume module



Condensation module

# Consumables



D2 lamp D2-100 (T60U)



D2 lamp D2-200 (T70, T80, T90, T92)



W lamp W-100 (T60V, T60U)  
W lamp W-200 (T70, T80, T90, T92)



Glass and Quartz cuvettes

**We offer a complete range of cuvettes**

# A3



The A3 Atomic Absorption Spectrometer is a high performance automated instrument designed to meet the requirements of the modern laboratory. Due to its versatility and performance it can be used for a wide range of applications including:

- Agricultural
- Food
- Geological
- Clinical
- Metal
- Petrochemical
- Environmental
- Mining
- Pharmaceutical

The versatile instrument is available in three configurations:

**A3F** – The instrument is equipped with a flame atomiser only. The positioning of which is fully controlled by the embedded computer system and AA Win software.

Three flame options are available to the user with the Air/Acetylene being the standard configuration. This flame can be used for nearly all standard elements while the N<sub>2</sub>O/Acetylene and the Air/LPG (Natural gas) are available as an option for the more demanding of elements. All three flame configurations offer coded burner for full safety protection.



**A3G** – The instrument is equipped with a Graphite Furnace Atomiser only. The graphite head is fixed into the optical path to maximise performance and eliminate drift.

The temperature of the transversely heated graphite tube is accurately controlled by means of a precision feedback system. Pyrolytically coated platform tubes are supplied as standard to improve the performance and eliminate many analytical problems associated with this technique.

**AA3FG** – The instrument is equipped with both Flame Atomiser and Graphite Atomiser as described above. Both configurations are installed into the instrument and can be changed over by a simple selection in the versatile AA Win software.

## FEATURES & FUNCTIONS

- Embedded PC system built into the instrument as standard on all configurations. Pre-installed AAWin3 software, user manuals, cook book and Windows operating system.
- AAWin3 software provides full control of the instrument and autosampler with easy method change for each technique. **Now GLP compliant.**
- Automatic 8 Hollow Cathode lamp turret controlled and optimised by the AAWin3 software. Operating lamp current and warm-up lamp current can be individually controlled to eliminate drift commonly associated with lamp warming.
- D2 lamp background correction system fitted as standard to all configurations. High energy D2 lamp and adjustable beam splitter mirror are optimised by the AAWin3 software.
- Self Reversal background correction system fitted as standard to all configurations. The high performance background system uses the same hollow cathode lamp as installed for the analysis. Minimum extra components are required and optical alignment is very simple. Self Reversal can be used for any element at any wavelength making it extremely versatile.
- High precision minimal optics ensures maximum light throughput to the computer controlled Czerny-Turner monochromator.
- A universal autosampler is available as an optional accessory which is conveniently mounted on the front of the A3 instrument.

- Optional XYZ autosampler provides high sample volume for multi element analysis along with scalable sample capacity for high ample throughput.
- Absorption and Emission modes are standard features in the AAWin3 software as well as peak height, peak area, sequential and manual integration modes.

## FLAME ATOMISER FEATURES

The flame atomiser offers three flame options:

Air/acetylene is the standard configuration with the N<sub>2</sub>O/acetylene and Air/LPG as options.

### Air/Acetylene

- The Air/Acetylene flame uses a 100mm single slot burner for the standard configuration.
- The high sensitivity (Cu 2ppm >0.280abs) is due to the efficiency of the fixed position glass nebuliser fitted as standard. An acid resistant replacement is available as an option.
- The flame can be easily set from blue lean flame through stoichiometric to fuel rich by means of computer control.

### N<sub>2</sub>O/Acetylene

- The N<sub>2</sub>O/Acetylene flame uses a 50mm single slot burner and is available as an optional extra.
- This flame configuration is used to measure elements less prone to ionization such as: Aluminium, Tin, Titanium, Calcium, Vanadium and Molybdenum.
- Switching from Air/Acetylene to N<sub>2</sub>O/Acetylene to Flame Off is fully controlled by the AAWinPro software.





# A3 continued



## Air/Propane (LPG)

- This flame uses a 3 slot burner and with the low pressure requirement it is also much safer to operate.
- Due to the low temperature of the flame it is ideal for analysing alkali metals such as: Potassium, Sodium and Lithium, especially when used in the emission mode.
- Some remote areas of the world have difficulty obtaining Acetylene gas of a high enough purity to operate the flame correctly, LPG can give a real alternative and offer comparable results for most elements throughout the wavelength range.

## Safety Features

- Pressure monitoring for all gases
- Burner Identification
- Flame sensor
- Drain Trap level Sensor
- Gas Leak Detector
- Over Pressure in Premix
- Safety Cut off Switch

## GRAPHITE ATOMISER FEATURES

The integrated Graphite Furnace Atomiser is available in two instrument configurations.

- In the A3G instrument the graphite furnace head is fixed into the light path so alignment with the optical path is simple and accurate.
- In the A3FG instrument the graphite furnace head is fixed behind the flame atomiser assembly and is motorised into position by a simple operation in the

AAWinPro software. The positions for the flame and graphite are saved making it easy to swap between modes for different analysis.

- The temperature of the transversely heated graphite tube is accurately controlled by means of a precision feedback system and has been designed to reduce analytical problems normally associated with this type of technique.
- Pyrolytically coated graphite tubes are used as standard and are manufactured to improve performance as well as increase the analytical life.
- Platform graphite tubes are supplied as standard and will accept volumes up to 20µl. Non-platform graphite tubes are also available as an optional extra.
- Up to 10 heat stages are available for the programming of the graphite atomiser. These can be set and stored within the AAWinPro software.
- The graphite tube is held in position by means of a gas piston. Replacement of the graphite tube is performed by a simple command in the AAWin Software.
- The graphite tube is efficiently cooled by an additional water circulation system (supplied separately).

## Safety Features

- Argon Gas pressure Sensor
- Water flow sensor
- Over Temperature Sensor
- Broken graphite tube protection Both configurations are installed into the instrument and can be changed over by a simple selection in the versatile AA-Win software.



# Specifications

Wavelength Range	185nm-910nm
Monochromator	Czerny-Turner Configuration
Spectral Bandwidth	0.1nm, 0.2nm, 0.4nm, 1.0nm, 2.0nm, (software selectable).
Wavelength Accuracy	± 0.15nm
Wavelength Reproducibility	< 0.05nm
Resolution	0.2nm ± 0.02nm
Baseline Stability	0.005A/30min
Sensitivity (Cu)	2 µg/ml Absorption > 0.28A (flame)
Detection Limit	Cu < 0.004 µg/ml (flame) Cd < 0.4 x 10 <sup>-12</sup> g (graphite furnace)
Repeatability	Cu < 0.7% (Air/Acetylene flame) Ba < 1.0% (Nitrous oxide/Acetylene flame) Cu < 2.0% Cd < 2.0% (Graphite Furnace)
Background Correction	Deuterium Arc, Self reversal
Characteristic Concentration	Cu < 0.02 µg/ml, Ba < 0.15 mg/ml (N2O/Acetylene)
Burner Heads	Titanium Alloy
Nebuliser	Inert Adjustable Nebuliser, Flow Rate from 2ml-6ml/min. Pt/Ir capillary available for high acid concentration
Atomization Chamber	Corrosion-resistant material
Position Adjustment	Automatic changeover (A3F) Manual (A3F) Automatic Setting of Optimum Height for Flame Burner.
Safety Functions	Burner Identification, Flame Sensor, Gas leak Sensor, Low Gas Pressure Sensor, Drain Trap Sensor, Power Loss Protection, Circulation Water (graphite), Over Temperature Sensor (graphite)

We reserve the right to modify, revise/upgrade, suspend or discontinue any Product in whole or in part, either temporarily or permanently, with or without notice.



# PF7



**The automated PF7 is an ideal instrument for elemental analysis in various market and research industries including:**

- Agricultural
- Pharmaceutical
- Food and beverage
- Geological
- Public health
- Metallurgical
- Clinical
- Petrochemical



## Standard Features

### Light source

- High intensity hollow cathode lamps for improved sensitivity and stability.
- Pre-aligned lamp assemblies for trouble free installation.
- All lamps are uniquely data coded offering important information to the PF Win operating software.
- Up to 3 lamps can be installed for simultaneous analysis.



### Optical System

- Double beam optical system to eliminate drift from the light source and the detector.
- Shielded optical design greatly reducing light interference.
- Enhanced signal to noise ratio for increased analytical sensitivity.
- Unique optical configuration for increased Fluorescence intensity. Up to twice the intensity found in traditional AFS systems.
- High Quantum Solar Blind detector fitted as standard to ensure optimum stability.

### Atomiser System

- High precision quartz tube designed for optimum performance, durability and long life.
- Adjustable height control for improved optimisation.
- Integrated 2 stage, fully sealed, fume exhaust system to decontaminate toxic elements and pollution.
- Gold Mesh fitted to the chimney removes any mercury pollutant.

### Hydride Generator

- Integrated continuous flow Hydride System.
- Gas pressure sampling offers maintenance free operation.
- Online auto dilution and multiple auto purge by gas driven sequential injection system.
- Fully sealed reservoir bottles for extended solution life.
- New design Gas Liquid Separator with magnetic stirring for improved repeatability of analytical results.
- Liquid Separator cooled directly by specially designed Peltier system to remove unwanted water in the formed hydride and greatly reduce Fluorescence quenching thus increasing the sensitivity.
- Unique high volume reagent storage positioned outside of the instrument to reduce contamination.
- Connection of carrier and reducer liquids to instrument using long life chemical resistant FEP tubing.

### Electronic Control

- High technology electronics and PCB components.
- PF Win 3.0 software offers full control of PF7 instrument and accessories.
- Windows operating software
- New features include: QC functions, online data sharing, self diagnostics, result and resource management.
- Full GLP version available for multiuser group management and log.





# PF7 continued

## Detection Limits

Element	Detection Limit (ug/l)	RSD %
Arsenic (As)	<0.01	<1%
Bismuth (Bi)	<0.01	<1%
Cadmium (Cd)	<0.001	<1%
Germanium (Ge)	<0.05	<1%
Mercury Hg)	<0.001	<1%
Lead (Pb)	<0.01	<1%
Antimony (Sb)	<0.01	<1%
Selenium (Se)	<0.01	<1%
Tin (Sn)	<0.01	<1%
Tellurium (Te)	<0.01	<1%
Zinc (Zn)	<1.0	<1%



Auto Sampler

## Accessories

### Auto Sampler

- X, Y, Z drive configuration.
- Fully controlled by PFWin software.
- Inert robust probe and FEP tubing.
- Improved probe wash. Simultaneous inner and outer wall wash.
- Large volume standard stock solution.
- Removable inert sample tray and rack.
- 3 sizes of sample racks available to accommodate 10ml, 25ml and 50ml test tubes.

### Speciation Analyser

- Built in HPLC pump. Isocratic/gradient
- Fitted with conventional column (optional column oven available).
- Manual or automatic sample input ( with optional auto sampler).
- Detect and separate inorganic and organic compounds.
- High separation performance
- Fast analysis < 12min.

## Speciation Detection Limits

Element	Specification	Detection Limit (ug/l)	RSD %
As	Arsenite (As III)	0.04	<5%
	Dimethylarsenic acid (DMA)	0.08	<5%
	Monomethylarsenic acid (MMA)	0.08	<5%
	Arcenate (As V)	0.2	<5%
Hg	Inorganic (Hg II)	0.05	<5%
	Methylmercury ( MeHg)	0.05	<5%
	Ethylmercury (EtHg)	0.05	<5%
	Phenylmercury (PhHg)	0.1	<5%
Se	Selenocysteine (SeCys)	0.3	<5%
	Selenite (Se IV)	0.1	<5%
	Selnomethionine (SeMet)	2.0	<5%
	Selenate (SE VI)	0.5	<5%
Sb	Sb III	0.1	<5%
	Sb V	0.5	<5%

# Specifications



Specifications	PF7
<b>Sample Atomisation</b>	
Atomiser	Quartz furnace tube with auto ignition
Furnace Heating	Computer controlled heating
Hydride Generator	Continuous flow high performance for cold vapour Mercury and hydride determinations of As, Se, Te, Bi, Sb, Sn, Zn, Pb, Cd, Ge
Gas Requirement	High purity Argon gas (99.99% ) 30psi
Exhaust System	2 stage filtration to decontaminate pollutants
<b>Sample and delivery</b>	
Carrier and Reagent Delivery	Gas pressure driven system
Gas/Liquid Separator	High efficiency gas/ liquid separator with magnetic stirring and Peltier cooling
<b>Optics</b>	
Optical Design	Short focal length non- dispersive double beam
Light Source	3 channel simultaneous element analysis using computer controlled modulated and pulsed hollow cathode source
Baseline Stability	< 1.5%
Baseline Noise	< 1.5%
Linear Range	>10 <sup>3</sup>
<b>Operating System</b>	
Interface	USB, RS232
Operating Software	PFWin software
Safety	Gas alarm for low pressure and flow
Power Requirements	100 – 240V 50Hz/60Hz 300VA
weight	Net: 35Kg Gross 50Kg
Dimensions	60cm x 57cm x 45cm





**Persee Analytics, Inc.**

Represented by  
Giangarlo Scientific Company, Inc.  
162 Steuben Street  
Pittsburgh, Pa. 15220

Web: [www.giangularscientific.com](http://www.giangularscientific.com)

Phone: 412-922-8850

Email: [sales@giangularscientific.com](mailto:sales@giangularscientific.com)

